


CANDO® Finger / Shoulder Ladder


- Offers progressive mobility of arm at shoulder, elbow and wrist joints
- 36 steps at 1.25" (3.2 cm) intervals; lowest step is 1.375" (3.49 cm) above the bottom edge of the ladder
- Step area: 2" (5.1 cm) wide x 44.75" (113.7 cm) long by 1.375" (3.5 cm) thick
- Total area: 2" (5.1 cm) wide x 47" (119.4 cm) long by 1.375" (3.5 cm) thick



Mounting Instructions

- Mount ladder to wall surface using two (2) fasteners (hardware not included) through the two (2) predrilled holes.
- It is suggested that the bottom edge of ladder is mounted 22" (55.9 cm) above the floor. At this suggested distance the top of the ladder is 69" (175.3 cm) from the floor. Check for adequate clearance before installing. This distance can vary depending on the desired range of motion.
- Make sure unit is securely fastened to the wall before putting it into service.
- FEI does not provide mounting hardware nor hardware recommendations due to unknown wall and building construction.

Please consult an installation professional who is familiar with your installation requirements to ensure compliance with applicable codes.



FEI Fabrication Enterprises Inc.
250 Clearbrook Rd, Suite 240
Elmsford, NY 10523 (USA)
tel: +1-914-345-9300 • 800-431-2830
fax: +1-914-345-9800 • 800-634-5370
FabEnt.com



EC REP

AJW Technology Consulting GmbH
Breite Strasse 3
40213 Düsseldorf (Germany)