



EXERCISE BAND RACK, HORIZONTAL, 5-ROLL



No	Item	Quantity
1	Bottom Casing Board	1
2	Separating Board	4
3	Round Acrylic tubing (15mm*2mm)	20
4	Round POM Rod (10mm)	4
5	Back Mounting Strip	1
6	Side Casing board	2
7	Wooden Peg (5mm)	2
8	Countersunk Screw (6mm*25 mm)	13
9	Top Casing Board	1
10	Plastic Cap (10mm*8mm)	13
11	Expansion Bolt (6mm)	4
12	Open End Wrench (#10)	1
13	Plastic Anchor Expansion Bolt (4mm*30mm)	4



1. Insert the four separating boards (2) onto slots of the bottom casing board (1), press them to enable all locate well and align on the same plane.
2. Insert four round POM rods (4) and twenty round acrylic tubing (3) onto holes of the four separating boards (2) properly as to the drawing (The cutting corners of four separating boards (2) are at top position.)
3. Assemble the four POM rods (4) with the two side casing boards (6).
4. Fasten the two side casing boards (6) with the bottom casing board (1) by four countersunk screws (8).
5. Put the back mounting strip (5) to sit on the cutting corners of four separating boards (2) and fasten it with the two side casing boards (6) by two countersunk screws (8).
6. Fasten the top casing board (9) with the back mounting strip (5) by three countersunk screws, and then fasten it with the two side casing boards (6) by four countersunk screws (8).
7. The plastic caps (10) are to cover all countersunk screw heads.
8. The four expansion bolts (11) are used for fastening the casing onto side wall by #10 open end wrench (12)
9. The four plastic anchor expansion bolts (13) are used to mount the casing onto roof, if needed.