

Cable Tensiometer









Intended Purpose

The **Baseline®** cable tensiometer is perfect for muscle strength measurement, job task analysis and functional capacity evaluation (FCE). It can be used to measure strength of various muscles.

Each unit includes two chains. One chain is attached to the unit and secured to a stationary object. The other chain is attached to the opposite side of the unit and to a pull cuff that is then pulled to measure strength. Maximum strength readings are held until unit is reset. Strength can be viewed in pounds or kilograms.

Each model includes two 5' chains and four pull cuffs that can attached to the ankle, thigh, wrist, and head.

Available in 3 models:

 12-0410
 Adult
 (660lb / 300 kg force capacity)

 12-0411
 Adolescent
 (330lb / 150 kg force capacity)

 12-0412
 Child
 (165lb / 75 kg force capacity)

MARNING:

- · Verify platform stability before use
- · Check chain and handle security
- Ensure proper patient positioning
- · Do not exceed device capacity
- · Clean all contact surfaces
- · Check all connecting parts

Breite Strasse 3

· Clean grip surfaces







40213 Dusseldorf Germany



NL-IM-00000248
Prinses Margrietplantsoen 33, Suite 123
2595 AM. The Haque (NL)

Tolerances:

The 660 lb cable tensiometer (12-0410) has an accuracy of ±10 lbs through its entire range.

The 330 lb cable tensiometer (12-0411) has an accuracy of ±8 lbs through its entire range.

The 165 lb cable tensiometer (12-0412) has an accuracy of ±5 lbs through its entire range.



12-0429 Head harness



Pull Cuffs

12-0410

Adult, cable

tensiometer unit

Wrap pull cuff around the body part being tested and secure using the hook and loop fastener. Attach pull cuffs to the Baseline® cable tensiometer chain using a snap or threaded oval.

12-0420 Ankle **12-0427** Shoulder

12-0422 Thigh **12-0428** Foot Stirrup, Large

12-0424 Wrist

12-0429 Head Harness

Disposal Method:

Dispose of item in accordance with the local/regional/national/international regulations.