

- non-slip material solves many stabilization and gripping problems encountered by people with limited dexterity, restricted mobility and reduced wrist or upper limb strength
- grips on both sides to prevent movement
- ideal for providing support, stability and confidence during everyday activities and exercise



color - Specify a color by adding the color abbreviation from the list at the end of each style, for example: 50-1500R is a red original 8"x10yd roll. **IF NO COLOR IS SPECIFIED, BLUE WILL BE DELIVERED.**

original Dycem® on rolls



- effective anti-slip material
- original material is 0.4mm with a thick non-woven center
- can be cut to any size and is malleable so it can be fitted around various objects
- *Y-yellow, R-red, B-blue, BLK-black, W-white, G-forest green*

50-1500	8" x 10 yard roll
50-1501	8" x 6½ foot roll
50-1502	8" x 3¼ foot roll
50-1508	16" x 16 yard roll
50-1505	16" x 10 yard roll
50-1506	16" x 6½ foot roll
50-1507	16" x 3¼ foot roll

grip-it Dycem® on rolls



- 0.6mm thick polyester-film base material
- more rigid than the original and should be used where a stiffer surface is required, such as under writing paper
- *F-floral, GH-gingham, C-clear*

50-1510	8" x 10 yard roll
50-1511	8" x 6½ foot roll
50-1512	8" x 3¼ foot roll
50-1515	16" x 10 yard roll
50-1516	16" x 6½ foot roll
50-1517	16" x 3¼ foot roll

comfort Dycem® on rolls



show with optional roll rack dispenser

- 1.1mm thick with center carrier
- soft foam texture has been designed to give extra padding during a variety of everyday activities and can help make therapy more relaxed and comfortable
- *blue only; shown with optional dispenser*

50-1521	8" x 6½ foot roll
50-1522	8" x 3¼ foot roll

roll rack dispenser

- fits up to 16"x 16 yards of any Dycem®

50-1509 roll rack dispenser

self-adhesive discs



- put on any piece of equipment to provide a secure base
- prevents movement especially on angled surfaces
- simply peel off the backing from the disc and press on to object to be secured
- *black only*

50-1580	5/8" diameter, 12 each
50-1581	7/8" diameter, 12 each

self-adhesive strips




- **original:** Dycem® strips are made with an adhesive backing to permanently secure the strip; *B-blue, Y-yellow*

50-1560 16" x 1-½" strip (3 ea.)

- **comfort:** Dycem® self-adhesive strips are available with foam for enhanced cushioning, protection, and vibration absorption; *blue only*

50-1570 16" x 2¼" comfort strips (3 each)

self-adhesive material panels



- **original:** self-adhesive material permanently secures the Dycem® piece in place; *B-blue, W-white, R-red, Y-yellow*

50-1530 16" x 1 yard original SA roll

- **comfort:** self-adhesive material with foam for enhanced cushioning, protection, and vibration absorption; *blue only*

50-1540 8" x 1 yard comfort SA roll