

# I-BRESIS™

## POWER TO DELIVER

### ACCELERATING DRUG DELIVERY THROUGH HYBRID IONTOPHORESIS TECHNOLOGY

## THE RIGHT MEDICATION, THE RIGHT APPLICATION

- Medications with positive or negative polarity
- Flexible, time-saving treatment options
- Three treatment modes: I-Bresis, Standard, or Patch-Only

### SHORTER TREATMENT OPTIONS:

- The I-Bresis Mode powers a 3-minute Skin Conductivity Enhancement to help break down the skin's resistance, resulting in a 40-80 mA-min treatment in approximately 1-2 hours.

### CONVENIENT:

- The I-Bresis Patch provides fill-and-go efficiency. Built-in batteries plus dose controller activation allows quick and effective application.



### I-BRESIS HYBRID IONTOPHORESIS SYSTEM

PART NUMBER	DESCRIPTION
500060	I-Bresis Patch, Box of 6
1360	I-Bresis Charging Station
1361	I-Bresis Dose Controller




## I-BRESIS CHARGING STATION

The I-Bresis Charging Station keeps the dose controllers charged between uses. Simply remove a wireless dose controller whenever a patch requires activation.



## I-BRESIS DOSE CONTROLLER

The dose controller activates I-Bresis and Standard Modes. Delivery flexibility allows clinicians to select the appropriate treatment combination for each individual patient.



## I-BRESIS PATCH

Built-in batteries and ergonomic design allow fast application options and fill-and-go efficiency. The I-Bresis Patch can be used with or without the dose controller.

## ACTIVE DELIVERY TO ACCELERATE RECOVERY:

Applying direct current drives penetration by opening pathways to the skin<sup>1</sup>.

### ORDERING INFORMATION

IMAGE	PRODUCT #	PRODUCT NAME	"PKG / ORDER UOM"
	1360	I-Bresis Charging Station	1 EA / EA
	1361	I-Bresis Controller	1 EA / EA
	500060	I-Bresis Patch, Box of 6	6 EA / BX
	5000008	Electrodes, Disposable, Transq, LG	12 EA / BX
	5000009	Electrodes, Disposable, Transq, MD	12 EA / BX
	5000033	Electrodes, Disposable, Transq, SM	12 EA / BX
	5000044	Iontophoresis System, Companion 80	6 EA / BX
	5000050	Electrodes, Optima, Disposable, SM	12 EA / BX
	5000051	Electrodes, Optima, Disposable, MD	12 EA / BX
	5000052	Electrodes, Optima, Disposable, LG	12 EA / BX

This pricing is valid until 7/1/2020 unless otherwise specified. Exact freight and tax will be added at the time of invoice.