Another NEW product from the makers of Relief Pak[®] cold & hot packs



- Instantly activated
- Reusable
- Economical
- Safe when wet
- Long-lasting heat









Recharge!

Once the metal disk in a Hot Button[™] compress is snapped, the liquid in the pad crystallizes and becomes hot right before your eyes! The Hot Button[™] compress generates heat for 30 minutes to 2 hours depending on the pad's size and the outside temperature. To recharge the compress, boil it for approximately 10 minutes. The activation and recharging of a pad can be repeated at least 100 times.

Size

3¹/₂" x 5¹/₂" 15" x 5¹/₂"

17" x 11"

11" x 10"

Contact Linda for Your Great Pricing!

| ltem | Description |
|------------|-------------------|
| A. 11-1025 | Small Compress |
| B. 11-1026 | Standard Compress |
| С. 11-1027 | Neck Compress |
| D. 11-1028 | Oversize Compress |







FABRICATION ENTERPRISES INC. PO BOX 1500 WHITE PLAINS, NEW YORK 10602 (USA) TEL: 1-914-345-9300 FAX: 1-914-345-9800 Linda@Fab-Ent.com



Suggested List

\$ 7.50

\$12.50

\$17.50 \$15.00

AUTHORIZED CE REPRESENTATIVE RMS UK LTD. 28 TRINITY ROAD NAILSEA, SOMESET BS48 4NU (UK)



