

PHC51G



29.5" x 23.38" x 19.88" (H x W x D)

Single chamber compact warming cabinet with glass door, digital thermostat, and lock

Highlights:

Compact size makes it ideal for use on top of counters

Digital controls conveniently located for easy adjustment and temperature viewing

Factory installed lock allows you to restrict access to staff

Product Features:

Compact size	Small footprint is ideal for countertop placement
Digital thermostat	Electronic controls are conveniently located for easy temperature adjustment
Electronic readout	Large LED temperature display ensures your staff can easily monitor the interior temperature
120 to 150°F Temperature Range	Designed and set for hospital blankets
Double pane tempered glass door	Heat-safe glass door offers a full view of stored contents
Self-closing door	Door is equipped with self-closing mechanism and magnetic gaskets to prevent heat loss
Factory installed lock	Keyed lock for a secure interior
Hammered aluminum interior	Ensures easy cleaning inside
Flush back	Space saving design with easier cleanability

PHC51G Specifications:

Overview	
Nominal Height	30.0" (76 cm)
Height of Cabinet	29.5" (75 cm)
Nominal Width	24.0" (61 cm)
Width	23.38" (59 cm)
Nominal Depth	20.0" (51 cm)
Depth	19.88" (50 cm)
Depth with door at 90°	42.0" (107 cm)
Capacity	2.0 cu.ft. (57 L)
Door	Glass
Cabinet	Gray
US Electrical Safety	ETL
Canadian Electrical Safety	ETL-C
Amps	5.5
Voltage/Frequency	115 V AC/60 Hz
Weight	75.0 lbs. (34 kg)
Parts & Labor Warranty	1 Year
UPC	761101050379
Warming Cabinet	
Door Swing	RHD
Reversible	No
Shelf Type	Wire
Shelf Qty	1
Thermostat Type	Digital
Interior Height	18.0" (46 cm)
Interior Width	18.75" (48 cm)
Interior Depth	15.38" (39 cm)
Temperature Range	120 to 150°F

