

PTHC65GUCLHD



34.25" x 23.5" x 23.5" (H x W x D)

Undercounter single chamber passive heating blanket warmer with glass door, stainless steel interior and exterior, digital thermostat, and lock

Highlights:

- Designed and assembled in the USA
- No internal fans or moving parts ensure silent, safe, and efficient operation
- Stainless steel interior for durable use and easy cleaning

Product Features:

Undercounter Warming Cabinet	24" wide blanket warmer designed and constructed in the USA
90 to 140°F temperature range	Designed and set to keep blankets warm in hospitals, medical and physical therapy offices, spas, and other facilities
Passive heating	Fixed heaters located at the bottom of the cabinet with no fans or moving parts to minimize lint buildup for safer, more efficient use
Built-In capable	Front-breathing design lets you make the best use of space by installing your appliance under the counter
Full capacity	4.57 cu.ft. interior in a slim-fitting 24" footprint
Digital thermostat	Electronic controls are conveniently located externally for easy temperature adjustment
Electronic readout	Large LED temperature display ensures your staff can easily monitor the interior temperature
Silent operation	No moving parts to ensure a quiet performance
Double pane tempered glass door	Heat-safe glass door offers a full view of stored contents, with a stainless steel door trim to highlight the interior

PTHC65GUCLHD Specifications:

Overview	
Height of Cabinet	34.25" (87 cm)
Width	23.5" (60 cm)
Depth	23.5" (60 cm)
Depth with Handle	25.38" (64 cm)
Depth with door at 90°	45.75" (116 cm)
Capacity	4.57 cu.ft. (129 L)
Door	Glass
Cabinet	Stainless Steel
US Electrical Safety	ETL
Canadian Electrical Safety	ETL-C
Amps	5.5
Voltage/Frequency	115 V AC/60 Hz
Weight	140.0 lbs. (64 kg)
Shipping Weight	125.0 lbs. (57 kg)
Parts & Labor Warranty	1 Year
UPC	761101108063
Warming Cabinet	
Door Swing	LHD
Reversible	Factory Reversible
Adjustable Shelves	Yes
Shelf Type	Chrome
Shelf Qty	3
Thermostat Type	Digital
Interior Height	26.0" (66 cm)
Interior Width	19.38" (49 cm)
Interior Depth	17.38" (44 cm)
Level Legs	4
Dolly	Optional
Temperature Range	90 to 140°F

