



TheraTouch® SW1

Portable Shockwave



The TheraTouch SW1 is Richmar's first radial pressure wave device. This device uses radial pressure waves to deliver a high energy mechanical force to the body's tissue to treat various musculoskeletal conditions. The radial pressure waves are transmitted through an ergonomic handpiece and applicator that can effectively penetrate 35mm into the tissue. Shockwave treatment can help tissue activation, bone healing, reduce muscle aches and pains, increase blood circulation and more.

The TheraTouch SW1 has an intuitive and simple touch screen interface. It offers a wide variety of preset protocols that were carefully programmed according to ISMST (International Society for Medical Shockwave Treatment) Guidelines and clinical testing. In addition to the preset protocols, there is also a manual mode in which users can set their own parameters, choose specific applicators and save custom treatments to a favorites list. The SW1 comes equipped with everything needed to provide effective treatment; including multiple applicators that can be easily switched out depending on the needs of the patient, a sleek carrying case and a rechargeable battery with a three hour battery life, allowing for easy portability and treatment on the go.

FEATURES:

- 11 Preset Body Part Protocols – Wrist, Elbow, Knee, Shoulder, Neck, Thigh, Ankle, Upper Back, Lower Back, Lower Leg and Foot
- Favorites List
- Three Applicator Sizes
 - 15mm, 20mm, 35mm
 - Steel-Round Applicators
 - Compound-Round Applicators (Patients Who Are More Sensitive to Treatment)
- Ability to Measure Energy in Both MJ and Bars
- Trilingual Interface
- Three-Year Warranty

TECHNICAL SPECIFICATIONS:

- **Power Supply:** 100-240V, 50/60Hz, 200VA
- **Lithium-Ion Battery:** DC 22.2V, 2200mAh
- **Frequency:** 1Hz – 22 Hz
- **Energy:** 60 – 185mJ, 1-4 bars
- **Main Unit Dimensions:** 11.46" x 6.69" x 4.78"
- **Handpiece Dimensions:** 10.2" x 4.48"
- **Main Unit Weight:** 4.63 lbs

STANDARD ACCESSORIES:

SW1000-HP	SW1 Handpiece	1
SW1000-HPS	SW1 Handpiece Stand	1
SR15-APP	Steel Round 15mm Applicator	1
CR15-APP	Compound Round 15mm Applicator	1
CR20-APP	Compound Round 20mm Applicator	1
SWAPP-RT	Applicator Removal Tool	1
SWORING95	Shockwave O-Ring 9.5mm	10
SWORING15	Shockwave O-Ring 15mm	10
SWORING20	Shockwave O-Ring 20mm	10
SWORING35	Shockwave O-Ring 35mm	10
DQ8000X	Power Cord	1
SW1BATT	Battery Module	1
USB	User Manual	1

OPTIONAL ACCESSORIES:

SR20-APP	Steel Round 20mm Applicator
SR35-APP	Steel Round 35mm Applicator
CR35-APP	Compound Round 35mm Applicator

ORDER DETAILS:

SW1000	TheraTouch SW1
---------------	----------------

WARNING: This product can expose you to chemicals including Lead and Di(2-ethylhexyl)phthalate (DEHP), which are known to the State of California to cause cancer or birth defects or other reproductive harm. For more information, visit www.P65Warnings.ca.gov.

