



JetBoots Prime

Wireless.

Ultra-Portable.

Safe, Science-Backed
Compression Therapy.

Recharge your legs anytime, anywhere with the ultra-portable JetBoots Prime compression boots, which weigh only 6 lbs and fold up compactly to fit in a backpack. Featuring a fully wireless design with no wires or hoses, JetBoots Prime offer hassle-free recovery wherever you are. Utilizing TruGrade Technology™, these boots accelerate recovery, reduce muscle soreness, and improve performance by delivering precise, science-backed pressure levels sequentially up the leg. The advanced built-in pump system ensures consistent and balanced treatment, optimizing circulation and effectively removing metabolic waste more than competitors.

Science-Backed Benefits

- Accelerate recovery
- Reduce muscle soreness
- Improve performance
- Improve flexibility & mobility
- Boost circulation

Key Features

- Fully wireless and lightweight- no wires, no cables, no cords
- TruGrade Technology™ delivers safe precise pressure levels for optimal circulation
- Simple one-touch control panel and LED indicators built right into the boots for easy treatment customization.
- 4 adjustable pressure levels and 4 time settings for customized treatments
 - Pressure levels: 25, 50, 75, 100 mmHg
 - Time settings: 20 min, 40 min, 60 min, and continuous
- Made with easy-to-clean medical-grade materials to withstand heavy use and prevent bacteria build-up
- Up to 3 hours of recovery with ability to use while charging
- Available in 3 sizes: Short, Regular & Long
- Includes drawstring travel bag
- Lighter and more portable, only 6 lbs.



 Therabody

PERFORMANCE

\$499
USD

\$599
CAD