



JetBoots PRO Plus

Wireless.

LED Light +
Vibration Therapy.

Safe, Science-Backed
Compression Levels.

The most advanced recovery boots on the market, JetBoots PRO Plus, use cutting-edge wireless pneumatic compression, vibration therapy, and infrared LED light to enhance your recovery. These features work together to improve circulation, strengthen muscles, relieve pain, reduce inflammation, and speed up healing. You can either customize your sessions or use one of the 8 expert-designed programs for optimal results. JetBoots PRO Plus are engineered to boost performance and accelerate recovery, giving you a competitive edge.

Scientifically-Proven Benefits

- Accelerate recovery
- Improve athletic performance
- Relieve muscle & joint pain
- Maximize warm-up
- Boost healing
- Increase strength

Key Features

- Powerful combination of pneumatic compression, vibration and infrared LED light therapies for faster, more effective recovery
- Wireless design with built-in pumps for maximum convenience
- FastFlush Technology delivers compression cycles faster than competitors with full inflation and deflation in as little as 60 seconds.
- TruGrade Technology delivers safe, precise compression sequentially up the leg through 4 overlapping pressure chambers to optimize circulation
- 8 preset treatments designed by leading sport scientists
 - Recovery, Warm-Up, Lower Leg, Upper Leg, Knee, Joint Therapy, Pain Relief, and Sleep Prep
- Built-in LCD screen with one-touch controls for quick and easy operation
 - Customizable time intervals: 10 – 60 min in increments of 5 min
 - Customizable pressure levels: 20 – 100 mmHg in increments of 5 mmHg
 - 3 vibration intensity levels: low, medium & high
 - LED: on / off
- Hygienic, non-porous material and interior lined design with one smooth, easy-to-clean surface
- Battery life ranges from 2.5- 4 hours (varies by therapies and settings utilized)
- USB-C charging and ability to use while charging
- Available in 3 sizes: Short, Regular & Long



IS Pneumatic
Compression
Therapy



Wireless
Design



LED Light
Therapy



Vibration
Therapy

 Therabody

PERFORMANCE

\$999
USD