

CaviCide1™ & CaviWipes1™

Surface Disinfectants

CaviCide1 ready-to-use liquids and CaviWipes1 disinfecting wipes are fast and effective with a 1-minute contact time and 1-step cleaning/disinfection for all pathogens listed.*



Features & Benefits



1-Minute Kill Time

Decreased risks of cross-contamination



1-Step* Cleaning and Disinfecting

Allows for faster room turnover



Clinically Relevant Efficacy Claims

Kills organisms that pose a threat in today's environment



Conveniently Packaged

Ready to use, provides flexible infection prevention at your fingertips



Materials Compatibility

Compatible with common hard, non-porous surfaces found in the clinical environment

Ordering Information

CaviCide1 Ready-to-Use Liquid

Part No.	Description
13-5024	24 ounce spray
13-5000	1 gallon bottle
13-5025	2.5 gallon bottle
13-5055	55 gallon drum

CaviWipes1 Pre-Saturated Wipes

Part No.	Description
13-5100	6" x 6.75", 160-count
13-5150	XL 9" x 12", 65-count
13-5155	Singles 9" x 12", box of 50
13-5224	Flat Packs 7" x 9", 45-count
13-1175	Canister Wall Bracket

*Always follow Product Label for Cleaning/Disinfection Directions.

1-Minute Efficacy Against:

Bacteria

*Acinetobacter baumannii**

Staphylococcus aureus

Pseudomonas aeruginosa

Salmonella enterica

Klebsiella pneumoniae

Enterobacter cloacae

Burkholderia cepacia

Fungi

Candida albicans

Candida auris

*Trichophyton interdigitale**

Drug Resistant Bacteria

Multi-drug resistant (MDR)

Acinetobacter baumannii

Methicillin Resistant

Staphylococcus aureus (MRSA)

Methicillin Resistant

Staphylococcus epidermidis (MRSE)

Vancomycin Resistant

Enterococcus faecalis (VRE)

Vancomycin Intermediate

Staphylococcus aureus (VISA)

ESBL *Escherichia coli*

Carbapenem-Resistant *Klebsiella pneumoniae* (CRKP) (CRE)

*CaviCide1 only

CaviCide1 EPA Reg. No. 46781-12

CaviWipes1 EPA Reg. No. 46781-13

Viruses

Human Immunodeficiency Virus (HIV-1)

Hepatitis B Virus (HBV)

Hepatitis C Virus (HCV)

Herpes Simplex Virus Type 1 & 2

Human Coronavirus

Influenza A

SARS-CoV-2

Measles virus

Norovirus*

Rotavirus*

Mycobacterium

Mycobacterium tuberculosis (TB)



Scan QR Code to learn more
or visit us at [metrex.com](https://www.metrex.com)