

Inject-Ease[®]

Injections Made Easy

Conquer the anxiety of injections with the Inject-Ease[®], the device that makes injections easy. Simply place your loaded syringe in the Inject-Ease, place the tip against your skin, and press the button to automatically deliver the needle through the skin. You control the rate at which the medicine is injected. The special tip is designed to help reduce the pain of injection. Making injections easier helps reduce anxiety - an independent clinical study showed that 91.6% of participants reported less fear of injections with the Inject-Ease¹.

From the people at AMBIMEDINC, makers of world-class medical device products.



Key Features/Benefits

- Complete needle coverage - Reduces fear of injection by hiding the needle
- Special tip - Designed to help reduce pain through pressure anesthesia
- One-button ease of use - Provides independence to those with limited dexterity and coordination
- Designed for one-handed use - Allows access to harder to reach injection sites
- One-press needle insertion - Insert the needle through the skin automatically and easily
- Designed for BD* insulin syringes - works with 30, 50 or 100 unit syringes
- Comes with spacer rings for customization of injection depth
- Convenient for travel - maintains preloaded syringe integrity

Key Selling Points

- Industry-leading 1 year warranty
- Usable with subcutaneous self-injection treatments, including diabetes and other conditions²
- AMBIMEDINC-recognized quality
- Affordable

¹ Clinical study available at www.ambimedinc.com/injections/iestudy.htm

² Consult your physician regarding use of the Inject-Ease with your condition.

Please note that not all products are available in all markets. For full details, please contact your local reseller or AMBIMEDINC.



(actual size)

Specifications

Length: 5.4 in. (13.72 cm)

Width (at widest point): 1.7 in. (4.31 cm)

Injection Type: Subcutaneous

Syringe Compatibility:

Most insulin and allergy syringes with fixed needles are compatible with the Inject Ease. Some syringe examples are BD, Terumo, Monoject, Ulticare, Easy Touch, Relion and Sure Comfort. The Inject-Ease is also compatible with some prefilled syringes with medication such as Humira by Abbvie. The 1ml BD 5/8 in. Tuberculin (TB) syringe is also compatible when a large adapter-spacer ring is used (must order adapter-spacer ring separately).

- 100 unit (1.0 ml, 1cc) BD* insulin syringe;
- 50 unit (0.5 ml, 1/2cc) BD* insulin syringe; or
- 30 unit (0.3 ml, 1/3cc) BD* insulin syringe.

Package includes: 2 spacer rings – 1/8 in. (0.317 cm) and 1/4 in. (0.635 cm) – to decrease injection depth.

Warranty:

AMBIMEDINC warrants that each new Inject-Ease shall be free of all defects and guarantees normal performance for a period of 1 year.

Available Packaging Configurations



Inject-Ease is a trademark of AMBIMEDINC. All rights reserved.

*Other trademarks are property of their respective holders.

PN/REF: IE400

Contents: (1) Inject-Ease, 3 adapter tips (30, 50, & 100 unit), 1/8 in. (0.317 cm) and 1/4 in. (0.635 cm) spacer rings, clear plunger cap, and instructions for use.

UPC Code: 0-25294-87400-7

AMBIMEDINC
World Class Healthcare Products

Specifications are subject to change without notice.
© 2022 AMBIMEDINC. All rights reserved.