



# Manosplint Digital Heat Pan

- Perfect thermoplastic heating solution for hospital or clinic settings
- High-quality stainless-steel unit features a bi-fold lid that easily slides open and closed with a built-in insulated handle
- For safety, a low water sensor automatically stops the pan from heating if the water level is too low
- Removable perforated steel tray keeps materials in place
- Four clear LED indicators include:
  - heating indicator
  - over-temperature alarm
  - timer indicator
  - set temperature indicator
- Heat pan is thermostatically controlled and includes a timer with an automatic start function, so you can preheat the pan before you begin fabricating
- **9 Liter model** is an efficiently sized pan that allows for portability while its smaller volume means there is a quicker heating time for the water
- **28 Liter model** is a standard sized pan that is deep enough for all day fabricating
- **56 Liter model** is a larger sized pan that accommodates larger patterns and sheets. It is just as deep as the 28 Liter model, but longer and wider

24-6545 Manosplint Digital Heat Pan, 9 Liter  
24-6546 Manosplint Digital Heat Pan, 28 Liter  
24-6547 Manosplint Digital Heat Pan, 56 Liter

