

CORE PRODUCTS® TROCHANTER BELT

- Provides maximum stability and support to the hips and pelvis
- Helps relieve pain and improve hip joint stability
- Easily adjust tension to ensure comfort and support
- Durable hook and loop closure for adjustability and convenience
- Provides relief for various hip, pelvis and lower back conditions
- Customized fit for enhanced comfort



CUSTOMIZED FIT

Available in one size to fit most users, this belt is a versatile option to accommodate a variety of body shapes and sizes. Customize the belt by simply cutting excess fabric to the ideal length for the individual.

Trim excess fabric using a sharp scissors. Be sure to leave a minimum of 6 inches of material beyond the metal loop for leverage and sizing accommodations.

Available From:



PROVIDES SUPPORT TO TROCHANTERIC REGION

Trochanter Belts are primarily used to address pain and discomfort within the trochanteric region. The trochanteric region consists of the greater and lesser trochanter (top of thighbone or femur), as well as various muscles, tendons, and ligaments within the hips and pelvis. By reducing excessive movement and stabilizing the hips, this belt provides support to the area and helps reduce strain to the surrounding tissues. Help correct pelvic tilt or rotation by using this belt to encourage proper posture and correct alignment that may put stress on the joints, muscles, and tissues in the hip region.

The Trochanter Belt is designed to comfortably stabilize the hips and pelvis to help ease pain. This belt helps to limit excessive movement or rotation to provide relief for various hip, pelvis, and lower back conditions, such as: trochanteric bursitis, hip arthritis, hip dysplasia, and postpartum recovery. This belt can also be used to help improve posture and alignment, injury rehabilitation, or for support during high-impact or repetitive physical activities.

COMFORT AND STABILITY

Durable, plush, 2" wide hook and loop material provides maximum support to the hips and pelvis. This belt contours to the natural curves of the hips for comfort and stability. The low profile, flexible design moves with you during normal daily activities such as walking, running, or standing for extended periods. This Trochanter Belt allows the user to make a one-time adjustment for size; no need to fit each use.

SPECIFICATIONS

Product Selection: Core Products®
Trochanter Belt, Black LSB-6147

Product Weight: .18 lbs

Dimensions: 2" x 67"

Materials: Nylon, Neoprene

Case Quantity: 6

Case Dimensions: 12" x 12" x 6"

Case Weight: 2 lbs

Latex: Caution: This product contains natural rubber latex which may cause allergic reactions

Country of Origin: Assembled in the USA of globally sourced materials

California Residents Only:

WARNING: Cancer and Reproductive Harm – www.P65Warnings.ca.gov

*Stabilized
Motion* 