

## Always Available Active Floor Recovery

Despite our best efforts, falls remain a serious threat—especially in high-pressure environments like hospitals, assisted living communities, nursing homes, care facilities, and emergency scenes. In tight spaces or urgent situations, traditional lifting aids can be bulky, heavy to carry, slow to assemble, or simply unavailable.



The **FloorRiser Active™** is your fast, stable, and **solo-operable solution**. Built for rehab teams, acute care environments, and first responders, it's designed for safe floor recovery—anywhere, anytime.

### Solo Responder Operation

One-button inflation for quick, safe lifts—designed for clinicians and first responders in urgent, tight, or wet spaces. Made from medical-grade, durable, and easy-to-clean materials for acute care environments.



### No Assembly Needed

Ready to deploy in seconds—just position and inflate.



### Compact & Portable

Weighs just 14 lbs (6.5 kg) and comes with a lightweight backpack for easy transport.



### Powerful Performance

With a 600 lbs weight capacity, it lifts individuals up to 29 inches using 7 inflatable layers—delivering maximum stability and support when it matters most.



### Battery Powered

Equipped with 2 rechargeable batteries and a 2-year full coverage warranty. Features automatic deflation for added convenience.

