

Rediscover walking ease, comfort and freedom.



SAFER BY DESIGN



Patented UPWalker



Imitation Product

ENGINEERED FOR SAFETY

SAFETY & STABILITY BY DESIGN

The patented UPWalker is the only upright mobility device designed from the ground up to reduce lateral and forward tip-over fall risk. Tested and certified to ISO standards (ISO 11199-3: 2005) for indoor and outdoor use.



The New UPWalker Lite is a lightweight, compact mobility solution designed to enhance daily living.

- The UPWalker Lite is a smaller, lighter, lower-priced version of the original UPWalker, especially well-suited for indoor use.
- Designed to maneuver in tighter spaces and for users who might be challenged by a heavier walker.
- Weighing just 15.5 pounds, the UPWalker Lite is easy to lift and transport.
- Increases mobility and independence for the elderly and mobility challenged.
- Enables users to stand upright and look ahead with better support, greater dignity and confidence to remain more active.
- Clinical study of LifeWalker technology showed improved posture, safer, more stable feeling, and reduced pain in legs, back, hands and wrists.
- Includes a beverage holder, backrest and personal item bag.



UPWalker® Lite

* Limited Lifetime Warranty covers frame against defects in materials and workmanship.

UPWalker®, LifeWalker® and Protostar® are registered trademarks of Protostar, Inc.,
dba LifeWalker Mobility Products

The UPWalker® Lite is designed with simple features for ease of use.

1. Handgrip
2. Parking/hand brake
3. Armrest pad
4. Height adjustment tube
5. Height adjustment lock
6. Sit-to-stand handles
7. Rear wheel
8. Brake adjustment knob
9. Front wheel
10. Seat
11. Seat rail
12. Backrest
13. Personal item bag
14. Beverage holder



Push-button height adjustment.



Folds easily for transport or storage.



UPWalker Lite Product Specifications

Features	Specifications
User Height Range	4'7"-5'10" (140-178 cm)
User Weight Capacity	300 lb (136 kg)
Length	29" (74 cm)
Width	24.8" (63 cm)
Armrest Height Range	32.5"-43.5" (82.5-110.5 cm)
Seat Width & Depth	18" (43 cm) W x 10" (25 cm) D
Seat Height From Ground	21" (53 cm)
Sit-to-Stand Handle Height	27" (68.6 cm)
Unit Weight	15.5 lb (7 kg)
Wheel Size	8" Front; 6" Rear
Folded Dimensions	
Length	29" (74 cm)
Height	40" (101.5 cm)
Width	8.5" (21.5cm)
Warranty	Lifetime Warranty on frame; 6-months on non-durable components & accessories.

Safe patented design. For detailed patent information go to www.lifewalker.us/patents-and-trademarks.