



A healthier,  
safer and  
more  
dignified  
walking  
experience.

Upright Posture Walker designed to be safer & more comfortable than traditional walkers and rollators.

## SAFER BY DESIGN



Patented UPWalker



Imitation Product

ENGINEERED FOR SAFETY

## SAFETY & STABILITY BY DESIGN

The patented UPWalker is the only upright mobility device designed from the ground up to reduce lateral and forward tip-over fall risk. Tested and certified to ISO standards (ISO 11199-3: 2005) for indoor and outdoor use.



The revolutionary UPWalker® restores activity and independence for the elderly and mobility challenged.

- Enables users to stand upright and look ahead with better support, greater dignity and confidence to walk farther and remain more active.
- Preliminary clinical studies revealed improved walking posture, safer, more stable feeling, and reduced pain in legs, back, hands and wrists.
- Helpful in treatment and mobility for people with Orthopedic, Cardiovascular, Neurological, Pulmonary and other health disorders.
- Stylish, easy-folding design, lightweight for transporting and storage, simple armrest adjustment for the perfect height in seconds.
- Cushioned sit-to-stand handles ease transition to standing position.
- Standard accessories include drink holder, seat backrest, personal bag.
- Sturdy frame with user weight capacity up to 350 lbs (Large model).



\* Limited Lifetime Warranty covers frame against defects in materials and workmanship.

# UPWalker®

P.O. Box 721929, San Diego, CA 92172 • (877) 488-0822 • LifeWalkerMobility.com

UPWalker®, LifeWalker® and Protostar® are registered trademarks of Protostar, Inc.,  
dba LifeWalker Mobility Products



# UPWalker® features provide benefits not available in traditional walkers.



## Adjustable-Height Armrests

Padded armrests enable users of varying height to stand tall & look ahead with confidence and dignity.



## Lockable Brakes

Ergonomic handbrakes are easy to use for better control, stopping and parking.



## Comfortable Seat

Fabric seat with backrest support slides front to back for sitting or walking. (Naugahyde seat standard on Large units.)



## Sit-To-Stand Handles

Padded handles make transition between sitting and standing easier and safer.



## Multi-Terrain Wheels

Large 8-inch rubber wheels help to move smoothly over bumps and uneven surfaces.



## Sturdy, Foldable Frame

Stylish, lightweight, folds easily for transport and storage yet strong enough to support 350 lbs (Large units).

## UPWalker Accessories:

- **Backrest Support**  
Detachable backrest provides comfortable support while sitting.
- **Beverage Holder**  
Attaches to handgrip arm for easy access to your water bottle.
- **Personal Item Bag**  
Polyester bag with fold-over top to conveniently stow personal items.
- **Luxury Item Bag**  
Elegant design, sturdy material, vented side pockets for portable oxygen concentrators (POC) and easy-access zippered top.
- **Shopping Bag**  
Perfect for groceries; features handles to use as a carry bag.
- **Smartphone Holder**  
Adjustable to securely hold phones of various sizes for handsfree operation.
- **Cane/Umbrella Holder**  
Keep a cane or umbrella within easy reach with this two-piece holder.
- **Flashlight/Taillight**  
Serves as a headlight, flashing safety light and rear-facing taillight.

The UPWalker is available in 2 sizes to accommodate user height from 4'7" to 6'5" and weight up to 350 lbs (Large).

## UPWalker Product Specifications

FEATURES	STANDARD (H200)	LARGE (H200-L)
<b>Product Specifications</b>		
User Height Range	4'7"-5'10" (140-178 cm)	4'10"-6'5" (147-196 cm)
User Weight Capacity	300 lb (136 kg)	350 lb (159 kg)
Length	33.5" (85 cm)	36" (91.5 cm)
Width	23.5" (60 cm)	25.5" (65 cm)
Armrest Height Range	32.5"-43.5" (82.5-110.5 cm)	36"-49.5" (91-126 cm)
Seat Width & Depth	18" (43 cm) W x 11" (28 cm) D	20" (51 cm) W x 15" (38 cm) D
Seat Height From Ground	20" (51 cm)	21" (53.5 cm)
Sit-to-Stand Handle Height	30" (76 cm)	32" (81 cm)
Unit Weight	23.5 lb (10.5 kg)	25.5 lb (11.5 kg)
Wheel Size	8"	8"
<b>Folded Dimensions</b>		
Length	33.5" (85 cm)	36" (91.5 cm)
Height	37.5" (95 cm)	40" (101.5 cm)
Width	10.5" (27 cm)	10.5" (27 cm)
<b>Warranty</b>		
Lifetime Warranty on frame; 6-months on non-durable components & accessories.		

Safe patented design. For detailed patent information go to [www.lifewalker.us/patents-and-trademarks](http://www.lifewalker.us/patents-and-trademarks).