


AIR-GO™ WALKER



Walk on Air™

The Air-Go Walker redefines mobility—delivering comfort, ease, and freedom beyond traditional walkers. It's comfortable on the hands, stable underfoot, and ready when it's time to rest. With oversized, cushioned Airgonomic™ handles, water-dispersing treaded wheels, and a padded mesh seat, it gives users both freedom of movement and confidence in every step. For you, it delivers showroom differentiation, stronger customer loyalty, and reliable margins with MAP pricing you can count on.

The support of a walker meets the convenience of a rollator



Comfortable Hold, Every Time

Oversized Airgonomic™ handles cushion hands and are easy to grasp — even for users with arthritis or reduced dexterity.



Easily Rest

Enjoy the freedom to sit and rest when needed on the ¾-length open-and-go padded mesh seat—designed shorter to let users walk naturally—while four height adjustments ensure a comfortable, customized fit.



Compact, Stable, Travel-Ready

The scissor-fold frame collapses in seconds for storage and transport, and the 6" water-dispersing treaded wheels with rear ski glides provide smooth, stable movement across various terrain.



Available Accessories



Storage Bag
Item # 102662-BAG



Cup Holder
Item # 102662-CUP



Product Specs

Height Range	Adjustable for users 4'9" - 6'7"
Product Dimensions	25.6" (W) x 25.2" - 26.9" (L) x 33.1" - 36" (H)
Seat Height	20" - 23"
Caster Dimensions	6" (D) x 1" (W)
Folded Dimensions	7.5" (W) x 26.2" (L) x 20.5" (H)
Product Weight	8 lb
Weight Capacity	300 lb
Warranty	Limited Lifetime
Retail Box Dimensions	21.5" (L) x 7.1" (W) x 28.5" (H)

Item #	Description	UPC
MD30100BK	Air-Go Walker, Onyx	822383031071
MD30100BL	Air-Go Walker, Sapphire	822383031088

