

WORLD'S LIGHTEST WHEELCHAIR!

Feather™ CHAIR FOR KIDS



Just when we thought we couldn't make a chair any lighter! The Feather Chair for Kids weighs only 12 pounds and is the world's lightest pediatric wheelchair and is available in 12" and 14" wide. Our newest chair has a weight capacity of 200 pounds and has available quick-release wheels and anti-tippers.

At only 12 pounds – this chair is easy to lift when traveling with children. A great chair for getting out and about.

The Feather Kids Wheelchair comes in 12" and 14" wide seats and includes colorful, washable nylon fabric and a seat belt.

SPECIFICATIONS

Weight Capacity	200 lbs	Front Caster Size	6"
Product Weight (total)	17 lbs	Rear Wheel Size	20"
Product Weight (w/out Wheels)	12 lbs	Axle Type	Quick-Release
Seat Widths Available	12", 14"	Armrest Height (from Seat)	8"
Seat Depth	13"	Armrest Type	Fixed Height, Non-Removable
Seat-to-Floor Height	16"	Back Height	15"
Overall Width	21.5", 23.5"	Upholstery Type	Nylon
Overall Height	36"	Upholstery Color	12" (Orange), 14" (Blue)
Folded Depth	27"	Warranty on Frame	One Year Limited
Folded Height	25"	Warranty on Wearable Parts	Six Months Limited
Folded Width	12"		

Lightweight



Foldable



Locking Brakes



Arrives Assembled

