## GANDO 🕱 Balance Board Combo™

## Economical choice - Each board does the job of 3!







- · Each board comes with 2 instability balls
- Screw in 2 instability balls to create a rocker board, or screw in 1 instability ball to make a wobble or asymmetric wobble board
- Use to balance train and strengthen the lower body
- Durable, heavy-duty and non-skid, yet lightweight and compact for portability
- Circular board has 3 instability ball placement options
- Rectangular board has 5 instability ball placement options
- · Each board has handles/slots to add band or tubing

Circular Board (16" diameter) w/ 2 Instability Balls								
instability ball difficulty			tilt	each				
yellow	x-easy	1" height	7°	10-2020				
red	easy	1.5" height	11°	10-2021				
green	moderate	2" height	14°	10-2022				
blue	difficult	2.5" height	18°	10-2023				
black	x-difficult	3" height	21°	10-2024				

## Rectangular Board (15" x 18") w/ 2 Instability Balls

instability ball difficulty			tilt 15"	tilt 18"	each
yellow	x-easy	1" height	7°	7°	10-2030
red	easy	1.5" height	11°	10°	10-2031
green	moderate	2" height	15°	13°	10-2032
blue	difficult	2.5" height	18°	16°	10-2033
black	x-difficult	3" height	22°	19°	10-2034

## Additional balls and boards available

Instability Balls only									
		challenge	height	each					
yello	N	x-easy	1"	10-1760					
red		easy	1.5"	10-1761					
gree	n	moderate	2"	10-1762					
blue	)	hard	2.5"	10-1763					
blac	k	x-hard	3"	10-1764					
5 ball	set	10-1765							



instability balls

 
 Boards only

 10-1776
 16" diameter (hard)

 10-1786
 20" diameter (medium)

 10-1796
 30" diameter (easy)

 10-1756
 rectangular board (15"L x 18"W)





 FABRICATION ENTERPRISES INC

 PO Box 1500, White Plains, NY 10602 (USA)

 tel:
 800-431-2830 / +1-914-345-9300

 fax:
 800-634-5370 / +1-914-345-9800

 FabEnt.com
 info@fab-ent.com