

CANDO® Balance Pad

- Closed-cell foam provides excellent elasticity and energy absorption
- Coarse surface for traction on bottom and top surfaces
- Sanitized coating resists sweat and saliva and protects against microbes and odors
- Stack balance pads to increase balance difficulty
- Great for individual or group use
- Measures: 16"L x 20"W x 2.5"H

32-1500B pad (ea)

32-1500B-10 pad (10 ea)



FABRICATION ENTERPRISES INC
PO Box 1500, White Plains, NY 10602 (USA)
tel: 800-431-2830 / +1-914-345-9300
fax: 800-634-5370 / +1-914-345-9800
FabEnt.com info@fab-ent.com