



## Heating Units

- all units, including 4 & 6 pack versions, are **insulated**
- insulation prevents unit exterior from becoming too hot and causing burns
- heats moist heat packs in a water bath up to 170° F - use at 160° F for best results
- on/off switch illuminates green when unit is on
- heating light glows orange when unit is warming to desired temperature
- available with either 110V or 220V
- comes with removable plug set
- made from durable stainless steel for use in high-traffic clinics



4 Pack



6 Pack

Relief Pak® Heating Unit (110V)			Relief Pak® Heating Unit (220V)		
size	unit only	with packs	size	unit only	with packs
<b>Stationary</b>					
4 Pack	11-1981	11-1961	4 Pack	11-1991	11-1971
6 Pack	11-1982	11-1962	6 Pack	11-1992	11-1972
8 Pack	11-1983	11-1963	8 Pack	11-1993	11-1973
<b>Mobile</b>					
8 Pack	11-1984	11-1964	8 Pack	11-1994	11-1974
12 Pack	11-1985	11-1965	12 Pack	11-1995	11-1975
24 Pack	11-1986	11-1966	24 Pack	11-1996	11-1976

## Moist Heat Packs and Covers

- reusable
- use with Relief Pak® cover or towel
- maintains therapeutic moist heat for up to 30 minutes
- washable foam-filled terry cloth or all terry covers available for each size



Neck



Standard



Half Size

### Moist Heat Packs

	size	each	bulk (12 each)
Standard	10" x 12"	11-1310	11-1310-12
Neck	7" x 24"	11-1311	11-1311-12
Oversize	15" x 24"	11-1312	11-1312-12
Half size	5" x 12"	11-1313	11-1313-12
Spine small	10" x 18"	11-1314	11-1314-12
Spine large	10" x 24"	11-1315	11-1315-12
Knee-shoulder	10" x 20"	11-1316	11-1316-12
Circular	10" dia	11-1205	11-1205-12

### Foam-Filled Hot Pack Covers

Standard	20" x 24"	11-1360	11-1360-12
Neck	9" x 24"	11-1361	11-1361-12
Oversize	24" x 36"	11-1362	11-1362-12
Half size	12" x 15"	11-1363	11-1363-12
Circular	12" dia	11-1266	11-1266-12
Standard w/pocket	24" x 24"	11-1364	11-1364-12
Oversize w/pocket	24" x 36"	11-1365	11-1365-12

### All-Terry Covers

Standard	20" x 24"	11-1366	11-1366-12
Neck	9" x 24"	11-1367	11-1367-12
Oversize	24" x 36"	11-1368	11-1368-12
Half size	12" x 15"	11-1369	11-1369-12
Circular	12" dia	11-1227	11-1227-12

Visit [FabEnt.com/Relief-Pak](http://FabEnt.com/Relief-Pak) for our full line of hot and cold products