


BASELINE[®] EVALUATION INSTRUMENTS

Pneumatic Bulb Dynamometer with reset



2" squeeze bulb dynamometer measures PSI.

Red maximum force indicator holds maximum force reading until reset. Max indicator can be reset to zero by turning the center reset knob counterclockwise.

12-0291: 0-30 PSI, with max pointer, 0.5 PSI gradations

12-0293: 0-15 PSI, with max pointer, 0.5 PSI gradations

This product contains dry natural rubber.



Fabrication Enterprises Inc.
250 Clearbrook Rd, Suite 240
Elmsford, NY 10523 (USA)
tel: +1-914-345-9300 • 800-431-2830
fax: +1-914-345-9800 • 800-634-5370
FabEnt.com



AJW Technology Consulting GmbH
Breite Strasse 3
40213 Düsseldorf (Germany)