

REF 12-1030, 12-1031, 12-1032

BASELINE[®] **EVALUATION
INSTRUMENTS**

X-ray Goniometers 180 degree range



Intended Purpose:

Use to evaluate and document the joint range of motion and angles on X-ray negatives.

- 180° head has two opposing scales marked in 2½° increments
- Black anodized aluminum with easy-to-read white radiopaque markings

Available in 3 arm lengths:

12-1030 14" (35 cm)

12-1031 8" (20 cm)

12-1032 6" (15 cm)

⚠ WARNING:

- Clean and disinfect device before and after each patient use
- Verify measurement markings are clearly visible
- Ensure smooth operation of all moving parts
- Align device properly with anatomical landmarks
- Apply gentle pressure during measurements
- Stop if patient reports pain or discomfort
- Replace if accuracy becomes compromised
- Do not attempt to bend or modify the device

Measurement Tolerance

This device measures degrees with a tolerance of ±1 degree.

Disposal Method

Dispose of item in accordance with the local/regional/national/international regulations.

 Fabrication Enterprises Inc
250 Clearbrook Road, Suite 240
Elmsford, NY 10523 (USA)
tel: +1-914-345-9300 • 800-431-2830
fax: +1-914-345-9800 • 800-634-5370
FabEnt.com

 2862

EC REP

MD

AJW Technology Consulting GmbH
Breite Strasse 3
40213 Dusseldorf Germany

 **MedEnvoy**
NL-IM-000000248
Prinses Margrietplantsoen 33, Suite 123
2595 AM, The Hague (NL)

 Made in
United States

©2025 FEI, all rights reserved.

Baseline is a trademark of Goldberg. rev0325