




		★ GOOD ★	★★ BETTER ★★	★★★ BEST ★★★	★★★ BARIATRIC ★★★
Item #		P503	P703	P903	P750
Product Name		Prime Care Bed Model P503	Prime Care Bed Model P703	Prime Care Bed Model P903	Prime Care Bed Model P750
					
Height Range		7.8" - 30"	7" - 30"	7" - 30"	11" - 30"
Weight Capacity		450 lbs	550 lbs	600 lbs	750 lbs
Resident Weight Capacity		425 lbs	500 lbs	550 lbs	700 lbs
Bed weight		215 lbs	215 lbs	222 lbs	282 lbs
Deck Width		35"	35"	36" - 39" - 42"	39" - 42" - 48"
Width option		None	Opt Integrated 39", 42"	Integrated Push Button	Integrated Push Button
Length		76" - 80"	76" - 80" - 84"	80" - 84" - 88"	80" - 84" - 88"
Length option		--	Option 88"	STD Integrated	STD Integrated
Mobility Heights		Any	Any	Any	Any
One Touch Floor Lock		--	✓	✓	✓
Frame Deck Design:		Grid	Slat	Slat	Slat
Vertical Lift		✓	✓	✓	✓
Motors		3	4	4	4
Trend /Rev-Trend		--	✓ - only w/ staff	✓ - only w/ staff	✓ - only w/ staff
Staff Control		--	Optional	Optional	Optional
Auto Contour		✓	✓	✓	✓
Leg/Heel Lift		--	✓	✓	--
Comfort (Cardiac)Chair		--	✓ - only w/Staff	✓	✓
Floor Lock Casters		✓	✓	✓	✓
All Wheels Swivel		✓	✓	✓	✓
Back Lite Pendent		--	--	✓	✓
Under Bed Light		--	--	✓	✓
Battery Back up		--	Optional	Optional	Optional
IV Pole		--	✓	✓	✓
Warranty	Frame	10 years	15 years	15 years	15 years
	Welds	15 years	Lifetime	Lifetime	Lifetime
	Electrical	2 years	5 years	5 years	5 years

LTC Bed Accessories		
Assist Bars		
PLTCAB	PLTCAR	PLTCDR
Assist Bar - P903/P703/P503/P750	Assist Rail - P903/P703/P503/P750	Deluxe Assist Rail - P903/P703/P750
		

LTC Bed Accessories			
P703 Accessories			
P703WE	P703LE	P703SC04	P703-BBK
Width Extension 39"-42"	Length Extension - 88"	Staff Control	Battery Backup w/Bracket
			

LTC Bed Accessories			
Accessories			
P903WE	P903SC05	P903-BBK	P903TRP
Width Extension 39"-42" for P903/P750	Staff Control - Ti-Motion motors for P903/P750	Battery Backup w/ Bracket for P903/P750	Transporter for P903/P703/P503/P902/P601/P301
			

LTC Bed Accessories											
Bed Ends											
35BWBE-A	35BWBE-B	35BWBE-FM	42BWBE-A	42BWBE-B	42BWBE-FM	35MCBE-A	35MCBE-B	35MCBE-FM	42MCBE-A	42MCBE-B	42MCBE-FM
35" Brownlee Woods Bed End OSLO Cherry	35" Brownlee Woods Bed End OSLO Oak	P35" Brownlee Woods Bed End OSLO Figured Mahogany	42" Brownlee Woods Bed End OSLO Cherry	42" Brownlee Woods Bed End OSLO Oak	42" Brownlee Woods Bed End OSLO Figured Mahogany	35" Mill Creek Bed End OSLO Cherry	35" Mill Creek Bed End OSLO Oak	P35" Mill Creek Bed End OSLO Figured Mahogany	42" Mill Creek Bed End OSLO Cherry	42" Mill Creek Bed End OSLO Oak	42" Mill Creek Bed End OSLO Figured Mahogany
											