




Re-Order Number	Descriptions	Size	Units
553010	Bed Wedge - Cozy Cover	7.5"H x 27"D x 24"W	Pair
553011	Bed Wedge - Vinyl Only 15 Degrees	7.5"H x 27"D x 24"W	Each
553030	Bed Wedge - Cozy Cover	10"H x 27"D x 24"W	Pair
553031	Bed Wedge - Vinyl Only 20 Degrees	10"H x 27"D x 24"W	Each
553050	Bed Wedge - Cozy Cover	12"H x 27"D x 24"W	Pair
553051	Bed Wedge - Vinyl Only 25 Degrees	12"H x 27"D x 24"W	Each

Purpose

The multi-purpose foam wedges comfortably elevate your patient's head or feet. More effective than pillows or rolled towels, they resist slipping and sliding for secure positioning. Solid with or without covers, Skil-Foam Wedges are available in many sizes.

Cover Laundering Instructions

The cover may be laundered and dried at temperatures up to 180°F (82°C). Bleaching is not recommended. Bleach will, over time, destroy the thread and fibers and shorten the useful life of the cover.