


Low-Shear I cover reduces skin friction and shear. Fluid-proof.



Vinyl cover is incontinent-proof and wipes clean for easy care.

High-density foam base prevents bottoming-out and increases sitting comfort.

Re-Order #	Description	Size	Unit
751010	Gel-Foam Vinyl Cushion w/LSI Cover	16"W x 16"D x 2.5"H	EA
751011	Gel-Foam Vinyl Cushion w/LSI Cover	18"W x 16"D x 2.5"H	EA
751013	Gel-Foam Vinyl Cushion w/LSI Cover	18"W x 18"D x 2.5"H	EA
751012	Gel-Foam Vinyl Cushion w/LSI Cover	20"W x 16"D x 2.5"H	EA
751014	Gel-Foam Vinyl Cushion w/LSI Cover	20"W x 18"D x 2.5"H	EA
751018	Gel-Foam Vinyl Cushion Only	16"W x 16"D x 2.5"H	EA
751019	Gel-Foam Vinyl Cushion Only	18"W x 16"D x 2.5"H	EA
751021	Gel-Foam Vinyl Cushion Only	18"W x 18"D x 2.5"H	EA
751020	Gel-Foam Vinyl Cushion Only	20"W x 16"D x 2.5"H	EA
751039	Gel-Foam Thin-Line Vinyl Cushion w/LSI Cover	16"W x 16"D x 1.5"H	EA
751040	Gel-Foam Thin-Line Vinyl Cushion w/LSI Cover	18"W x 16"D x 1.5"H	EA
751044	Gel-Foam Thin-Line Vinyl Cushion w/LSI Cover	18"W x 18"D x 1.5"H	EA
751049	Gel-Foam Thin-Line Vinyl Cushion w/LSI Cover	20"W x 18"D x 1.5"H	EA
751042	Gel-Foam Thin-Line Vinyl Cushion w/LSI Cover	20"W x 16"D x 1.5"H	EA
914481	Gel-Foam Ultra-Thin Vinyl Cushion w/LSI Cover	16"W x 16"D x 1"H	EA
914482	Gel-Foam Ultra-Thin Vinyl Cushion w/LSI Cover	18"W x 16"D x 1"H	EA
914483	Gel-Foam Ultra-Thin Vinyl Cushion w/LSI Cover	20"W x 16"D x 1"H	EA
751025	Gel-Foam Vinyl Cushion w/Coccyx Cutout & LSI Cover	18"W x 16"D x 2.5"H	EA
754528	Gel-Foam Vinyl Cushion w/Coccyx Cutout & LSI Cover	18"W x 16"D x 3.5"H	EA
751038	Gel-Foam Thin-Line Vinyl Cushion Only	16"W x 16"D x 1.5"H	EA
751041	Gel-Foam Thin-Line Vinyl Cushion Only	18"W x 16"D x 1.5"H	EA
751045	Gel-Foam Thin-Line Vinyl Cushion Only	18"W x 18"D x 1.5"H	EA
751043	Gel-Foam Thin-Line Vinyl Cushion Only	20"W x 16"D x 1.5"H	EA
751200	Sittin Pretty Bimini Blue Gel-Foam Vinyl Cushion	16"W x 16"D x 2.5"H	EA
751201	Sittin Pretty Bimini Blue Gel-Foam Vinyl Cushion	18"W x 16"D x 2.5"H	EA
751202	Sittin Pretty Bimini Blue Gel-Foam Vinyl Cushion	18"W x 18"D x 2.5"H	EA



Gel-Foam Coccyx Cushion - 751025



Sittin Pretty Bimini Blue - 751200

Purpose

Skil-Care's Gel-Foam Cushion is a gel over foam designed to relieve pressure and increase comfort. The water based gel surface provides cooling action that reduces skin-damaging perspiration.

Easy to Use

Place cushion on a level surface and secure to the chair with safety ties.

Easy Care

The Low-Shear outer cover can be laundered and dried at temperatures up to 180°F (82°C). Do not use bleach.