


Footrest Bag

Reorder #	Unit	Wt.	Unit
□914362	14.5" W x 21"L	.75 lbs	Ea.



Purpose

Most wheelchairs have footrests that are removable and aren't interchangeable with other model wheelchairs. Often, when removed, the footrests are difficult to locate which necessitates a large time expenditure to locate. The Skil-Care footrest bag can be used on any wheelchair and avoids looking for the detached footrests.

Use Instructions

The vinyl footrest bags are to be placed on the back of the wheelchair by placing the loops around each of the respective wheelchair handles. The footrests should be placed in the bag bottomside first.

See Illustration.

Care Instructions

The vinyl footrest bags can be cleaned with a wipe clean detergent.