

Mettler® ultrasound

shared features

- small portable size
- hand held operation
- rechargeable battery and charger for portable power
- packages include rechargeable battery (13-3408), battery charger (13-3410), power cable (13-3411) and battery case (13-3409)



Sonicator® 718



Sonicator® 719

Sonicator® 718

- full clinical power, 2 W/cm² continuous, 3 W/cm² pulsed portable ultrasound
- 5.5 cm² soundhead

| | | |
|---------|------------------------|----------|
| 13-3400 | Sonicator 718® package | 2,265.00 |
| 13-3401 | Sonicator 719® package | 2,265.00 |

Sonicator® 719

- 3.3 MHz 0.9 cm² portable ultrasound

accessories

| | | |
|---------|---|--------|
| 13-3404 | protective cap for 718 applicator | 22.00 |
| 13-3405 | protective cap for 719 applicator | 22.00 |
| 13-3406 | AC adaptor 718/719 | 285.00 |
| 13-3407 | power cable for AC adaptor (7181), (110-120V) | 62.50 |
| 13-3408 | rechargeable battery | 155.00 |
| 13-3410 | battery charger | 155.00 |
| 13-3409 | battery case for 13-3407 | 57.50 |
| 13-3411 | power cable for battery charger | 14.50 |
| 13-3402 | unit stand for 718/719 | 32.50 |
| 13-3403 | pouch for 718/719 | 82.50 |

Mettler® neuro stimulation

Sys*Stim® 240

- 9 stimulation waveforms for versatility; optional laser or cluster light therapy
- waveforms may be combined with on-off programs for neuromuscular education
- frequency and amplitude modulation programs for pain management
- patented M-Wheel™ for easy navigation through menus
- complete anatomical library using built-in protocols with electrode placement
- add either the Laser (on page 129) or Cluster (infrared) probes for complete clinical versatility

| | | |
|----------|---|----------|
| 13-3430 | Sys*Stim® 240 | 3,180.00 |
| 13-3431 | laser and cluster applicator kit | 3,040.00 |
| 13-3433* | laser applicator kit 785nm, 80mW | 1,735.00 |
| 13-3435* | cluster applicator kit 600/950nm, 500mW | 1,400.00 |

* kit has applicator, holders, hardware and 2-glasses



Sys*Stim® 228

- very versatile unit
- 4 discreet waveforms including Interferential, Premodulation, Medium (Russian) and Biphasic
- additional versatility with various treatment modes: surge, reciprocation, surge with reciprocation and vector rotation
- use with any Mettler® ultrasound for combination therapy

| | | |
|---------|---------------|----------|
| 13-3440 | Sys*Stim® 228 | 2,000.00 |
|---------|---------------|----------|



Mettler® combination ultrasound and stimulation

Sonicator® Plus 920

- touch screen combination device offers Interferential (4-pole), Premodulated (2-pole), Medium frequency, Biphasic, Hi-Volt, TENS, Micro-current and DC waveforms
- allows clinician to treat 3 patients at the same time
- 61 pre-set protocols with additional space for your own protocols
- 11 languages and internal safety features
- additional .09 cm² applicator available

| | | |
|---------|---------------------|----------|
| 13-3450 | Sonicator® Plus 920 | 4,290.00 |
|---------|---------------------|----------|

