

patient lifts

Drive® hydraulic powered lifts

hydraulic with 4-pt cradle


- 450 lb weight capacity
- silver vein steel construction
- high-performance hydraulics raise or lower individuals gradually and safely
- standard with 5" casters
- 3" casters (2 locking and 2 non-locking) available for low bed use
- adjustable-width base




hydraulic with 6-pt cradle

- 450 lb weight capacity
- chrome-plated steel provides maximum strength
- high-performance hydraulics raise or lower individuals gradually and safely
- easy-to-operate caster brakes provide security
- includes sling chains



lift with 5" casters 
41-0100 with 5" casters 512.50
additional 3" casters
41-0101 3" casters (4 ea) 50.00


lift with 5" casters 
41-0110 with 5" casters 575.00
additional 3" casters
41-0111 3" casters (4 ea) 50.00

Drive® battery powered lifts

battery powered lift

- 450 lb weight capacity
- 6 point cradle
- fast, easy charge from AC outlet
- battery does not need to be removed to be charged
- for home-care use
- audible warning when battery power falls below 50%




lifts 
41-0120 patient lift 1,720.00
41-0121 w/rechargeable battery 1,795.00

accessories
41-0125 recharg. battery 390.00
41-0126 battery charger 407.50
41-0127 wall mount charger bracket 37.50

bariatric battery powered lift

- 600 lb weight capacity (can fit under bed)
- battery does not need to be removed to be charged
- for home-care use
- base legs adjust easily and lock securely into open position
- extra wide coated cradle with 360° rotation with 4-point hook up



lifts 
41-0130 patient lift 2,085.00
41-0131 w/rechargeable battery 2,620.00


accessories
41-0125 recharg. battery 390.00
41-0126 battery charger 407.50
41-0127 wall mount charger bracket 37.50

Lumex® patients lifts

sit-to-stand battery lift

- provides assistance to those who have trouble standing
- adjustable knee pad unit
- lower base height maximizes stability and fits under most beds
- warning beep for low battery charge and weight capacity levels




lifts 
41-0180 patient lift 2,660.00
41-0181 w/rechargeable battery 3,295.00

accessories
41-0175 recharg. battery 185.00
41-0176 battery charger 105.00
41-0177 wall mount charger bracket 105.00

everyday electric lift

- heavy-gauge steel construction
- six-point spreader bar with 360° rotation
- manually adjustable base
- 24V DC motor allows quiet smooth operation
- easy access emergency stop button
- 400 lb weight capacity




lifts 
41-0160 blue 790.00
41-0161 black/silver 790.00

easy lift battery patient lift

- heavy-gauge steel construction
- 400 lb weight capacity
- six-point spreader bar with 360° rotation
- lower base height maximizes stability and fits under most beds
- warning beep for low battery charge and weight capacity levels




lifts 
41-0170 patient lift 2,060.00
41-0171 w/rechargeable battery 2,060.00

accessories
41-0175 recharg. battery 185.00
41-0176 battery charger 105.00
41-0177 wall mount charger bracket 105.00

hydraulic 6-pt

- heavy gauge steel construction
- six point spreader bar with 360° rotation
- ergonomic placement of hydraulic pump handle facilitates ease of use
- disassembles for easy storage
- 400 lb weight capacity



lifts 
41-0150 blue 790.00
41-0151 black/silver 790.00
41-0152 black/silver with foot pedal 790.00