

MSS™ Tilt and Recline positioning system

- for children with moderate to severe physical involvement
- hydraulic actuator activates the tilt-in-space feature
- chair can be moved without having to reposition child

standard features:

- height, depth & angle adjustable multi-axis headrest
- adjustable seat depth: 3½"-12"D
- height & angle adjustable back: 14"-19"H
- adjustable length pelvic belt and H-harness
- height & width adjustable laterals, armrests, abductor and footrests
- removable & washable cushion covers
- 100 lb capacity; 35 lb product weight
- seat width 9"-10½"; overall width 18"; seat to footrest height 4"-10½"; seat to floor height 14"

MSS™ floor sitter

31-3710	MSS™ floor sitter	1,225.00
---------	-------------------	----------

MSS™ chair and mobile base

31-3711	MSS™ chair and low mobile base 17"H seat to floor	1,275.00
31-3706	low mobile base only	180.00
31-3721	MSS™ chair and high mobile base 24" seat to floor	1,500.00
31-3707	high mobile base only	180.00



MSS™ chair and low mobile base

Tumble Forms® corner chair

- large base prevents tipping; built-in abductor helps leg extension and reduces extensor thrust
- adjustable harness support and hip strap securely position child
- removable tray (22"W x 10"D x 10"H) can lie flat for feeding or tilt to 30°
- tray cutout 14"W x 8"D
- 60" maximum user height
- 125 lb capacity
- overall 22"W x 26"D x 20H; seat 13"W x 10"D x 20"H



30-3050 corner chair 687.50

Tumble Forms® square module seating system

- position children up to 60" height
- shapes modify the seating for a variety of settings
- 8-piece shape system contains 2 large rectangles, 2 small rectangles, 2 triangles, 1 ring and 1 disc
- 10-piece seating system includes 8 pieces above plus another disc and box seat



10-piece seating system



8-piece shape system (box seat not included)

30-3061 8-piece shape system 725.00
30-3062 10-piece seating system 897.50