

vestibular stimulation

plastic ball stool/trainer



plastic ball chair



CanDo® wobble ball



- versatile exercise ball chair holds a 55 cm ball (included) securely in place
- use as a chair or use as an inflatable exercise dome
- casters have locking mechanism to make chair stationary
- lock casters when using during vestibular exercises to prevent chair from rolling away
- ball color may vary

- mobile ergonomic seating system combines the benefits of ball therapy with the convenience of a chair base
- ball held securely in place
- removable back, without arms
- casters glide on carpet and hard surfaces
- includes ball; colors may vary
- light assembly required

- use for vestibular stimulation, core training and general balance exercises
- adjust the level of difficulty by adding or removing air from ball
- use the dome curved side down as a seat, rocker board, or to kneel on for vestibular activities
- with the flat side down, use for balance, core and strengthening exercises

30-1796 mobile exercise ball stool 100.00

30-1792	adult size w/ 55 cm ball	125.00
30-1795	child size w/ 35 cm ball	110.00
30-1797	locking casters, 2 each	15.00
30-1798	replacement ball-adult 55 cm	15.00
30-1789	replacement ball-child 35 cm	15.00

30-1941 wobble ball 60.00

CanDo® inflatable exercise balls

- use to improve balance, coordination, flexibility, strength and even just for fun
- used by therapists for vestibular stimulation
- non-slip surface is ribbed for extra security
- under inflate to give a soft, mushy feel or fully inflate for a firm, bouncy feel
- supports up to 300 lb (136kg)



	cm	in	color	
30-1800	30	11.8	blue	10.00
30-1801	45	17.7	yellow	12.50
30-1802	55	21.6	orange	15.00
30-1803	65	25.6	green	20.00
30-1804	75	29.5	red	25.00
30-1805	85	33.5	blue	30.00
30-1806	95	37.4	red	37.50

CanDo® inflatable saddle rolls



- easier to control than a ball because motion is limited to forward and backward direction
- aids in developing muscles, coordination and balance
- the "saddle" design offers extra support for the user
- supports up to 300 lb (136 kg)

	cm H x L	in H x L	color	
30-1725	40 x 90	16 x 35	yellow	30.00
30-1726	50 x 100	20 x 39	orange	40.00
30-1727	60 x 110	24 x 43	green	50.00
30-1728	70 x 120	28 x 47	red	60.00
30-1729	80 x 130	32 x 51	blue	70.00