


analog dynamometers for use with your existing software

300 lb hand dynamometer

100 lb pinch gauge

500 lb wrist dynamometer

500 lb push-pull dynamometer



- each dynamometer contains a transducer that generates an analog electronic signal; dynamometer does not produce digital signal
- the analog signal must be converted to a digital signal to be read by your computer software
- your in-house tech can do this with an A to D converter
- read information from your dynamometer dial face or directly on your computer with appropriate application software (not provided)
- computer compatible dynamometers are available in 2 styles:
 - 1) dial gauge with analog output
 - 2) analog output only

- Dual output: dial gauge with analog signal (see A)
- Single output: analog signal only (see B)

hand dynamometer

A	12-0021	300 lb (136 kg) dial gauge with analog output	775.00
B	12-0248	300 lb (136 kg) analog output only	575.00

pinch gauge

C	12-0023	100 lb (45 kg) dial gauge with analog output	700.00
	12-0238	100 lb (45 kg) analog output only	500.00

wrist dynamometer

D	12-0027	500 lb (226 kg) dial gauge with analog output	900.00
	12-0267	500 lb (226 kg) analog output only	700.00

push-pull dynamometer

E	12-0025	500 lb (226 kg) dial gauge with analog output	1,275.00
	12-0384	500 lb (226 kg) analog output only	1,075.00

Trade-in or upgrade your old dynamometer or pinch gauge

Now is the time to trade-in or upgrade your old dynamometer, pinch gauge or force gauge (any brand, any condition) for a new Baseline® measurement instrument.

trade-in program

trade in SKU	old product description	new product description	
12-0240T	any dynamometer	Baseline 200 lb hydraulic dynamometer	180.00
12-0235T	any pinch gauge	Baseline 50 lb hydraulic pinch gauge	120.00
12-0200T	any pinch gauge	Baseline 30 lb mechanical pinch gauge	110.00
12-0201T	any pinch gauge	Baseline 60 lb mechanical pinch gauge	110.00
12-0394T	any force gauge	Baseline 250 lb dial force gauge	412.50

upgrade program

trade in SKU	old product description	new product description	
12-0243T	any dynamometer	Baseline 200 lb hydraulic dynamometer, large dial	192.50
12-0246T	any dynamometer	Baseline 300 lb hydraulic dynamometer, ER™ large dial	207.50
12-0247T	any dynamometer	Baseline 300 lb hydraulic digital dynamometer	357.50
12-0239T	any pinch gauge	Baseline 50 lb hydraulic pinch gauge, large dial	152.50
12-0228T	any pinch gauge	Baseline 95 lb hydraulic pinch gauge, ER™ large dial	165.00
12-0237T	any pinch gauge	Baseline 100 lb hydraulic digital pinch gauge	357.50
12-0399T	any force gauge	Baseline 250 lb digital force gauge	495.00

Contact us for a return authorization number. A return authorization will be emailed to you. Include it with your old measurement instrument.