

Baseline™ push-pull dynamometers/manual muscle tester

- lightweight, small and portable
- easy to grasp with or without handles
- hydraulic, electronic or mechanical spring mechanism

muscle strength measurement

hand-held dynamometer lets you objectively measure manual muscle strength

job task analysis

measure actual push, pull and lift forces needed to perform a particular task (function)

functional capacity evaluation

quantitatively evaluate an individual's push, pull or lift capacity to perform a given task (function)

use for push or pull applications



muscle strength testing



job task evaluation



functional capacity evaluation

measure lifting



adjustable length chain

hydraulic

digital readout with push pad

analog readout with pull hook

1 year warranty

1 cm² circular pad



- comes with accessory pack: 3 push pads (padded curved, padded straight, and 1 cm² circular), 1 pull hook and 1 snap lock hook
- maximum remains until rest
- optional accessories available (see page 249)
- cushioned carrying case
- 1 year warranty; CE certified

electronic

with push pad

with pull hook

2 year warranty



- comes with accessory pack: 3 push pads, 1 pull hook, 1 snap lock hook, protective carrying case
- "load-cell" technology assures the highest standards of accuracy, repeatability and reliability
- maximum remains until reset
- optional accessories available (see page 249)
- 2 year warranty; CE certified

mechanical spring

with pull hook

with circular disc

with push pad



1 year warranty

- comes with accessory pack: 5 push pads: padded curved (x2), padded straight, and circular (x2); 1 pull hook and 1 snap lock hook
- maximum remains until rest
- 1 year warranty; CE certified

analog (dial) readout

12-0392	50 lb / 23 kg	600.00
12-0393	100 lb / 45 kg	675.00
12-0394	250 lb / 113 kg	750.00
12-0388	500 lb / 226 kg	825.00

digital (LCD) readout

12-0397	50 lb / 23 kg	900.00
12-0398	100 lb / 45 kg	975.00
12-0399	250 lb / 113 kg	1,050.00
12-0387	500 lb / 226 kg	1,125.00

12-0340	50 lb / 23 kg	900.00
12-0341	100 lb / 45 kg	1,000.00
12-0342	250 lb / 113 kg	1,100.00
12-0343	500 lb / 226 kg	1,200.00

12-0300	2.2 lb / 1 kg	150.00
12-0301	5.5 lb / 2.5 kg	150.00
12-0302	11 lb / 5 kg	150.00
12-0304	22 lb / 10 kg	150.00
12-0303	66 lb / 30 kg	150.00