


MicroFET 4™ grip/pinch tester



- combines grip and pinch strength testing in one device
- unit can perform four grip tests (maximum, 5-position, rapid exchange and fatigue) and three pinch tests (key, palmar and tip)
- 200 lb force capacity
- unit contains: test forms, wall chart, case, 3.7V rechargeable lithium ion battery and AC adapter

12-0298W	MF4, wireless	1,095.00
12-0298WC	MF4, with clinical software	1,595.00
12-0298WD	MF4, with data collection software	1,595.00
12-0298WCD	MF4, with both software packages	2,095.00

ErgoFET™ MMT



- displays peak force and time duration
- push attachments ensure proper muscle isolation for increased testing validity
- pull attachments
- push and pull forces accurate to within 1%
- 300 lb capacity
- includes: 3 push attachments, 2 pull hooks, D-ring strap, extension post, test forms and (2) 3.7V rechargeable lithium batteries with AC adapter

12-0461W	EF, wireless	1,395.00
12-0461WD	EF, with data collection software	1,895.00

MicroFET HandGRIP™



- ideal for routine screening and ongoing evaluations
- unit can perform four grip tests (maximum, 5-position, rapid exchange and fatigue)
- unit contains: test forms, wall chart, case, 3.7V rechargeable lithium ion battery and AC adapter

12-0277W	MF HandGRIP, wireless	355.00
12-0277WC	MF HandGRIP, with clinical software	855.00
12-0277WD	MF HandGRIP, with data collection software	855.00
12-0277WCD	MF HandGRIP, with both software packages	1,355.00

ErgoFET 500™ push-pull dynamometer



- measure peak force, average or 3-second average force
- push and pull attachments allow for a wide variety of strength measurement tests
- use removable handles when measuring higher forces
- 500 lb push capacity, 200 lb pull capacity
- includes: 3 push attachments, 2 pull hooks, D-ring strap, extension post, test forms, handles, wall chart, case, 3.7V rechargeable lithium battery and AC adapter

12-0460W	EF500, wireless	1,395.00
12-0460WD	EF500, with data collection software	1,895.00