


# AccuForce™ Preferred Colors®

# LATEX band

## AccuForce™ Preferred Colors® exercise band

- teach patients how much force to exert
- improve patient compliance at home
- same color sequence as Thera-Band® and CanDo® band

challenge level	yellow 1	red 2	green 3	blue 4	black 5
	easy		moderate		difficult



when a level of force is reached, the indicator shape will change from a rectangle to a square



50-yards



6-yards

### How much force is being exerted?

- each band has 4 rectangles
- when stretched, the rectangles progressively become squares
- read the force level directly from "square" on the stretched band

Percent of stretch	exerted force in lbs			
	50	100	150	200
1.00	1.75	2.50	2.75	
1.25	2.75	3.75	4.50	
2.25	3.50	4.50	5.25	
3.75	5.75	7.25	7.75	
4.50	7.00	8.50	9.75	



40 4-foot singles in dispenser box



4-foot singles

### 50-yard dispenser box

DISPENSER BOX	yellow	x-light	10-5921	85.00
	red	light	10-5922	92.50
	green	medium	10-5923	100.00
	blue	heavy	10-5924	107.50
	black	x-heavy	10-5925	115.00
	yel red grn blu blk	set (1 ea)	10-5928	475.00
		set w/ Dispens-a-Band® rack	10-5929	607.50

### 6-yard dispenser box

DISPENSER BOX	yellow	x-light	10-5911	12.00
	red	light	10-5912	13.00
	green	medium	10-5913	14.00
	blue	heavy	10-5914	15.00
	black	x-heavy	10-5915	16.00
	yel red grn blu blk	set (1 ea)	10-5918	65.00

### 40 4-foot singles in dispenser box

PRE-CUT SINGLES	yellow	x-light	10-5941	85.00
	red	light	10-5942	92.50
	green	medium	10-5943	100.00
	blue	heavy	10-5944	107.50
	black	x-heavy	10-5945	115.00
	yel red grn blu blk	set (1 ea)	10-5948	475.00

### 4-foot singles

PRE-CUT SINGLES	yellow	x-light	10-5901	2.75
	red	light	10-5902	3.00
	green	medium	10-5903	3.25
	blue	heavy	10-5904	3.50
	black	x-heavy	10-5905	3.75
	yel red grn blu blk	set (1 ea)	10-5908	15.25

## CanDo® accessory kit

includes accessories to enhance band or tubing workout

- standard foam handles - pair
- door disc anchor
- loop stirrup anchor
- use with all brands of band and tubing

10-5334	accessory kit	15.00
10-5330	foam handles with cinch, pair	10.00
10-5325	door disc anchor with cinch, each	5.00
10-5310	double loop stirrup anchor, each	1.75

### get more use from band and tubing with optional accessories



accessory kit



foam handle with cinch (pair)



double loop stirrup anchor



door disc anchor

Exercise band is made of natural rubber latex, which may cause allergic reactions • Thera-Band® colors licensed by Goldberg from Hygenic Corp. without endorsement

exercise