

bi-directional Endorphin™ units come standard with:

- pulse meter displays distance, time, speed, calories and pulse
- body stabilizing chest pad is adjustable to vertical or horizontal position with 7" range-of-motion

unit	bi-directional pedal	chest support	external power 110 volt	handles	programmable resistance	work levels WATTS	pulse meter	workout selection	screen size
e1	Y	Y	N	comfort	8 manual		Y	N	small
e2	Y	Y	Y	comfort	8 electronic plus manual	19-48	Y	Y	large
e3	Y	Y	Y	comfort	8 electronic plus manual	6-20	Y	Y	large
e4	Y	Y	Y	multi-position	8 electronic plus manual	19-48	Y	Y	large

UBE 300



- use on table or counter
- multi-position grip has 3 preset hand positions (22, 45 and 90°) to isolate and work different muscle groups
- 28"L x 16"W

comfort grip

[e1]	10-3600	945.00
[e2]	10-3605	1,595.00
[e3]	10-3610	1,595.00

multi-position grip

[e4]	10-3617	1,980.00
------	---------	----------

UBE 355 on stand



- adjustable wheelchair accessible stand
- crank arm adjustable from 36"-50" from floor allowing unit to be used in seated or standing position

comfort grip

[e1]	10-3620	1,815.00
[e2]	10-3625	2,425.00
[e3]	10-3630	2,425.00

multi-position grip

[e4]	10-3637	2,725.00
------	---------	----------

UBE 380 on stand with seat



- with height-adjustable platform
- semi-recumbent seat adjusts front to rear for proper user positioning; rotates 360° to provide easy access and egress
- seat assembly may be quickly removed for wheelchair accessibility

comfort grip

[e1]	10-3660	2,745.00
[e2]	10-3665	3,145.00
[e3]	10-3670	3,145.00

multi-position grip

[e4]	10-3677	3,495.00
------	---------	----------

LBE 300



- place unit on floor, in front of chair, bench or wheelchair
- 28"L x 16"W


foot pedal

[e1]	10-3640	1,045.00
[e2]	10-3645	1,695.00
[e3]	10-3650	1,695.00

foot plate

[e1]	10-3641	1,235.00
[e2]	10-3646	1,885.00
[e3]	10-3651	1,885.00

LBE + UBE 370



- comes with independent upper and lower ergometers
- height adjustable platform
- semi-recumbent adjustable seat

foot pedal/grip type

[e1]	10-3680	comfort grip	3,185.00
[e2]	10-3685	comfort grip	3,715.00
[e3]	10-3690	comfort grip	3,715.00
[e4]	10-3697	multi-position grip	3,995.00

foot plate/grip type

[e1]	10-3720	comfort grip	3,375.00
[e2]	10-3725	comfort grip	3,905.00
[e3]	10-3730	comfort grip	3,905.00
[e4]	10-3737	multi-position grip	4,185.00

accessories



accessories

10-2925	foot pedal	100.00
10-2921	foot plate	290.00
10-2919	comfort grip	25.00
10-2927	multi-position grip	550.00
10-2922	quad grip	80.00
10-2923	table	875.00