

Digi-Flex LiTE® Hand Exerciser

- Increase patient compliance at home
- Inexpensive option for patient giveaway or sale
- Individual finger locks
- Single-spring design
- Comes with exercise pamphlet
- Available in eight (8) resistance levels

force to compress
finger/hand (lb)

tan	xx-light	0.75	2.5
yellow	x-light	1.5	5.0
red	light	3.0	10.0
green	medium	5.0	16.0
blue	heavy	7.0	23.0
black	x-heavy	9.0	31.0
silver	xx-heavy	11.0	36.0
gold	xxx-heavy	13.0	43.0



[View Product Information Page](#)



Purchase Digi-Flex LiTE® from one of these suppliers. If your supplier can't full your order you can order directly from Fabrication Enterprises Retail Sales

Patterson Medical	Meyer PT	Mckesson Med-Surg	Scrip Hessco	FEI	Description	List	American 3B	Alimed	Direct Supply	Henry Schein	Medline
568971	FAB660TAN	943930	-	10-3770	Digi-Flex LiTE® - Tan (xx-light)	\$12.50	-	-	-	-	MDSP103770
568972	FAB660YEL	943931	-	10-3771	Digi-Flex LiTE® - Yellow (x-light)	\$12.50	-	-	-	-	MDSP103771
568973	FAB660RED	943932	-	10-3772	Digi-Flex LiTE® - Red (light)	\$12.50	-	-	-	-	MDSP103772
568974	FAB660GRN	943933	-	10-3773	Digi-Flex LiTE® - Green (medium)	\$12.50	-	-	-	-	MDSP103773
568975	FAB660BLU	943934	-	10-3774	Digi-Flex LiTE® - Blue (heavy)	\$12.50	-	-	-	-	MDSP103774
568976	FAB660BLK	943935	-	10-3775	Digi-Flex LiTE® - Black (x-heavy)	\$12.50	-	-	-	-	MDSP103775
568977	FAB660SIL	943936	-	10-3776	Digi-Flex LiTE® - Silver (xx-heavy)	\$12.50	-	-	-	-	-
568978	FAB660GLD	943937	-	10-3777	Digi-Flex LiTE® - Gold (xxx-heavy)	\$12.50	-	-	-	-	-
568981	-	-	-	10-3797	Digi-Flex LiTE® - Set of 5	\$57.50	-	-	-	-	MDSP103797
568982	FAB6605SET	-	-	10-3798	Digi-Flex LiTE® - Set of 5, with Plastic Stand	\$77.50	-	-	-	-	MDSP103798
568983	-	-	-	10-3799	Digi-Flex LiTE® - Set of 5 with Metal Stand	\$95.00	-	-	-	-	-
568979	-	943938	-	10-3778	Digi-Flex LiTE® - Set of 8	\$95.00	-	-	-	-	-
568980	FAB6608SET	943939	-	10-3779	Digi-Flex LiTE® - Set of 8, with 2 Metal Stands	\$170.00	-	-	-	-	-



Fabrication Enterprises Retail Sales PO Box 1076 White Plains, New York 10602
tel: 888-211-5242 / 914-407-0150
fax: 888-211-5765 / 914-407-0151
www.FEIRetail.com



AUTHORIZED CE REPRESENTATIVE
RMS UK LTD.
28 Trinity Road
Nailsea, Somerset BS48 4NU (UK)

Digi-Flex LiTE is a registered trademark of Goldberg © 2015 FEI all rights reserved