

Digi-Flex® Hand Exerciser

- The industry standard for 25 years
- For rehabilitation injuries or conditions affecting the hand and upper extremities
- Upper and lower spring design allows for more hand motion flexibility
- Develops isolated finger strength, flexibility, and coordination
- Available in eight (8) strength levels
- Comes with exercise pamphlet

force to compress finger/hand (lb)			
tan	xx-light	0.75	2.5
yellow	x-light	1.5	5.0
red	light	3.0	10.0
green	medium	5.0	16.0
blue	heavy	7.0	23.0
black	x-heavy	9.0	31.0
silver	xx-heavy	11.0	36.0
gold	xxx-heavy	13.0	43.0



[View Product Information Page](#)

Purchase Digi-Flex® from one of these suppliers. If your supplier can't full your order you can order directly from Fabrication Enterprises Retail Sales

Patterson Medical	Meyer PT	Mckesson Med-Surg	Scrip Hessco	FEI	Description	List	American 3B	Alimed	Direct Supply	Henry Schein	Medline
960617	FAB244TN	811082	-	10-0751	Digi-Flex® - Tan, xx-light	\$20.00	-	-	0-C8228	1222304	MDSP100751
A39715	FAB244YEL	766152	845013306	10-0740	Digi-Flex® - Yellow, x-light	\$20.00	W51124	5966/NA/15	0-27324	9603166	MDSP100740
A3973	FAB244RD	766153	845013304	10-0741	Digi-Flex® - Red, light	\$20.00	W51120	5966/NA/3	0-27325	9602203	MDSP100741
A3975	FAB244GRN	697425	845013303	10-0742	Digi-Flex® - Green, medium	\$20.00	W51121	5966/NA/5	0-37054	9608665	MDSP100742
A3977	FAB244BL	766154	845013302	10-0743	Digi-Flex® - Blue, heavy	\$20.00	W51122	5966/NA/7	0-37055	9603816	MDSP100743
A3979	FAB244BLK	771339	845013301	10-0744	Digi-Flex® - Black, x-heavy	\$20.00	W51123	5966/NA/9	0-37056	9607307	MDSP100744
960618	FAB244SIL	811083	-	10-0752	Digi-Flex® - Silver, xx-heavy	\$20.00	-	-	-	1222305	-
960619	FAB244GLD	811084	-	10-0753	Digi-Flex® - Gold, xxx-heavy	\$20.00	-	-	-	1222306	-
-	-	800284	845 0128	10-0745	Digi-Flex® - Set of 5	\$95.00	-	-	0-23173	9609271	MDSP100745
-	FAB242	539243	845 0129	10-0746	Digi-Flex® - Set of 5 with plastic rack	\$112.50	-	-	0-23174	9607610	MDSP100746
A39710	-	800285	-	10-0749	Digi-Flex® - Set of 5 with metal rack	\$132.50	-	-	0-23175	1155490	-



Fabrication Enterprises Retail Sales PO Box 1076 White Plains, New York 10602
 tel: 888-211-5242 / 914-407-0150
 fax: 888-211-5765 / 914-407-0151
 www.FEIRetail.com



AUTHORIZED CE REPRESENTATIVE
 RMS UK LTD.
 28 Trinity Road
 Nailsea, Somerset BS48 4NU (UK)

Digi-Flex is a registered trademark of Goldberg © 2015 FEI all rights reserved