

Digi-Flex Thumb® Hand Exerciser

- Thumbs can now be incorporated into the program
- Exercises fingers and thumb
- Design allows for both left and right thumb usage
- Thumb can attach in three different positions to accommodate all thumb sizes
- Thumb rotates on all planes
- Fingers and thumb buttons lock
- Comes with exercise pamphlet
- Available in eight (8) strength levels

force to compress
finger/hand (lb)

tan	xx-light	0.75	2.5
yellow	x-light	1.5	5.0
red	light	3.0	10.0
green	medium	5.0	16.0
blue	heavy	7.0	23.0
black	x-heavy	9.0	31.0
silver	xx-heavy	11.0	36.0
gold	xxx-heavy	13.0	43.0



[View Product Information Page](#)

Purchase Digi-Flex Thumb® from one of these suppliers. If your supplier can't full your order you can order directly from Fabrication Enterprises Retail Sales

Patterson Medical	Meyer PT	Mckesson Med-Surg	Scrip Hessco	FEI	Description	List	American 3B	Alimed	Direct Supply	Henry Schein	Medline
081655802	FAB661TAN	-	-	10-3760	Digi-Flex Thumb® - Tan (xx-light)	\$20.00	-	-	-	-	MDSP103760
081655810	FAB661YEL	-	-	10-3761	Digi-Flex Thumb® - Yellow (x-light)	\$20.00	-	-	-	-	MDSP103761
081655828	FAB661RED	-	-	10-3762	Digi-Flex Thumb® - Red (light)	\$20.00	-	-	-	-	MDSP103762
081655836	FAB661GRN	-	-	10-3763	Digi-Flex Thumb® - Green (medium)	\$20.00	-	-	-	-	MDSP103763
081655844	FAB661BLU	-	-	10-3764	Digi-Flex Thumb® - Blue (heavy)	\$20.00	-	-	-	-	MDSP103764
081655851	FAB661BLK	-	-	10-3765	Digi-Flex Thumb® - Black (x-heavy)	\$20.00	-	-	-	-	MDSP103765
081655869	FAB661SIL	-	-	10-3766	Digi-Flex Thumb® - Silver (xx-heavy)	\$20.00	-	-	-	-	-
081655877	FAB661GLD	-	-	10-3767	Digi-Flex Thumb® - Gold (xxx-heavy)	\$20.00	-	-	-	-	-
-	-	-	-	10-3769	Digi-Flex Thumb® - Set of 5	\$95.00	-	-	-	-	-
-	FAB6615SET	-	-	10-3785	Digi-Flex Thumb® - Set of 5 with plastic stand	\$112.50	-	-	-	-	-
-	-	-	-	10-3786	Digi-Flex Thumb® - Set of 5 with metal stand	\$132.50	-	-	-	-	-
-	-	-	-	10-3768	Digi-Flex Thumb® - Set of 8	\$150.00	-	-	-	-	-
-	FAB6618SET	-	-	10-3787	Digi-Flex Thumb® - Set of 8 with 2 metal stands	\$227.50	-	-	-	-	-



Fabrication Enterprises Retail Sales PO Box 1076 White Plains, New York 10602
tel: 888-211-5242 / 914-407-0150
fax: 888-211-5765 / 914-407-0151
www.FEIRetail.com



AUTHORIZED CE REPRESENTATIVE
RMS UK LTD.
28 Trinity Road
Nailsea, Somerset BS48 4NU (UK)

Digi-Flex Thumb is a registered trademark of Goldberg © 2015 FEI all rights reserved