



**FABRICATION
ENTERPRISES INC.**

To Whom It May Concern:

All items listed below can be cleaned/sanitized using soap and water, or a disinfectant wipe (such as Sanicloth, Super Sanicloth Wipes, Gym Wipes, CavieWipes, CaviCide1 Spray, Isopropyl Alcohol Wipes, Precision Optical Lens Cleaning Wipe, Quaternary Ammonia Wipes, ProSpray Wipes) or spray. Cleaner/sanitizer should be applied after every use. Bleach and any other commercial disinfectants, such as Oxivir, are not permitted. Using these can cause products to break down and will alter the quality:

- Airex Balance Pads
- Airex Exercise Mats
- Baseline Pulse Ox
- Baseline Finger Goniometer
- Baseline Tuning Fork
- Baseline Taylor/Babinski Reflex Hammer
- Baseline Hand Dynamometer
- Baseline Pinch Gauge (Mechanical & Hydraulic)
- Baseline Windmill Type Spirometer
- Baseline Dolorimeter
- Baseline Monofilaments
- CanDo Dumbbells
- CanDo Kettlebells
- CanDo Medicine Balls
- CanDo Low Powder Exercise Band
- CanDo Latex Free Exercise Band
- CanDo Foam Roller – White PE foam
- CanDo Foam Roller – Black Composite
- CanDo Foam Roller – Blue EVA Foam
- CanDo Exercise Balls (all versions)
- CanDo Peanut Rolls
- CanDo Twist-n-Bend
- CanDo Ball Stabilizer
- CanDo Wate Balls
- CanDo Massage Balls
- CanDo Band/Tubing Racks
- Classic/Deluxe Kickboards
- CanDo Hand Exercise Web Latex/Latex Free
- Cuff Weights
- CanDo Blue TufCoat®
- CanDo SoftGrip Hand Weights
- CanDo Dynamic Stretch Strap
- CanDo Donut Balls
- CanDo Inflatable Balance Disc
- Digi-Extend and Squeeze
- Digi-Flex
- Footstools
- Dellon Disk-Criminator
- Pouch Variable Wrist and Ankle Weight's
- Shoulder/Finger Ladder (plastic & wood)
- Skillbuilders Rolls
- Skillbuilders Wedges
- Transfer Boards – Wood/Plastic
- Hook and Loop material
 - o On products and stand alone
- Puttycise Theraputty Tools
- Rebounder - Adjustable, Circular, Square
- WaxWel Paraffin Bath

Sincerely,

Jason Drucker

SVP of Sales & Marketing